Cell Counting Kit-8

Cell Proliferation Assay and Cytotoxicity Assay

Technical Manual

for 100 tests, 500 tests, 1000 tests,
3000 tests, and 10000 tests

General Information

Cell Counting Kit-8 (CKK-8) allows very convenient assays by utilizing Dojindo's highly water-soluble tetrazolium salt. WST-8 [2-(2-methoxy-4-nitrophenyl)-3-(4-nitrophenyl)-5-(2,4-disulfophenyl)-2H-tetrazolium, monosodium salt] produces a water-soluble formazan dye upon reduction in the presence of an electron mediator, as shown in Figure 1.

CKK-8 is a one-bottle solution; no premixing of components is required. CKK-8, being nonradioactive, allows sensitive colorimetric assays for the determination of the number of viable cells in cell proliferation and cytotoxicity assays. WST-8 is reduced by dehydrogenases in cells to give an orange colored product (formazan), which is soluble in the tissue culture medium (Figure 2). The amount of the formazan dye generated by dehydrogenases in cells is directly proportional to the number of living cells.

Figure 3 shows that the cell proliferation assay using CKK-8 correlates well with the [³H]-thymidine incorporation assay. Thus, the CKK-8 assay can also be substituted for the [³H]-thymidine incorporation assay. As shown in Figure 4, the detection sensitivity using CKK-8 is higher than assays using other tetrazolium salts such as MTT, XTT, MTS or WST-1.


Figure 1. Structures of WST-8 and WST-8 formazan


Figure 2. Principle of the cell viability detection with Cell Counting Kit-8

Kit Contents

- 100 tests :	1 ml	x	1
- 500 tests :	5 ml	x	1
- 1000 tests :	5 ml	x	2
- 3000 tests :	5 ml	x	6
- 10000 tests :	100 ml	x	1

Storage

Store at 0-5°C

CKK-8 is stable over one year at 0-5°C with protection from light. Store it at -20°C for longer storage. Repeated thawing and freezing causes an increase in the background, which interferes with the assay. Please store the kit at 0-5°C for frequent use.

Required Equipment and Materials

- plate reader (450 nm filter)
- 96-well plate
- CO₂ incubator
- 10 µl and 100 - 200 µl multi-channel pipettes

General Protocol

Cell Number Determination

1. Inoculate cell suspension (100 µl/well) in a 96-well plate. Pre-incubate the plate in a humidified incubator (e.g., at 37 °C, 5% CO₂).
2. Add 10 µl of the CKK-8 solution to each well of the plate. Be careful not to introduce bubbles to the wells, since they interfere with the O.D. reading.
3. Incubate the plate for 1 - 4 hours in the incubator.
4. Measure the absorbance at 450 nm using a microplate reader.

To measure the absorbance later, add 10 µl of 1% w/v SDS or 0.1 M HCl to each well, cover the plate and store it with protection from light at room temperature. No absorbance change should be observed for 24 hours.


Figure 3 Correlation between CKK-8 assay and [³H]-thymidine incorporation assay.

Medium : HeLa.....MEM, 10% FBS
HL60.....RPMI1640, 10% FBS
Reagent : [³H]-Thymidine.....37 KBq/well
CKK-8.....10 µl/well
Incubation : [³H]-Thymidine assay..... 4 hours
CKK-8.....3 hours


Figure 4 Cell number determination using CKK-8 and other reagents.

Medium : HeLa.....MEM, 10% FBS
HL60.....RPMI1640, 10% FBS
Incubation : HeLa.....37°C, 5 % CO₂, 2 hours
HL60.....37°C, 5 % CO₂, 3 hours
Detection : CKK-8 (●) ... 450 nm, XTT (◆).....450 nm
MTS (▲).....490 nm, MTT (□).....570 nm

Cell Proliferation and Cytotoxicity Assay

1. Dispense 100 μ l of cell suspension (5000 cells/well) in a 96-well plate. Pre-incubate the plate for 24 hours in a humidified incubator (e.g., at 37°C, 5% CO₂).
2. Add 10 μ l of various concentrations of substances to be tested to the plate.
3. Incubate the plate for an appropriate length of time (e.g., 6, 12, 24 or 48 hours) in the incubator.
4. Add 10 μ l of CCK-8 solution to each well of the plate.
Be careful not to introduce bubbles to the wells, since they interfere with the O.D. reading.

5. Incubate the plate for 1 - 4 hours in the incubator.
6. Measure the absorbance at 450 nm using a microplate reader.

To measure the absorbance later, add 10 μ l of 1% w/v SDS or 0.1 M HCl to each well, cover the plate and store it with protection from light at room temperature. No absorbance change should be observed for 24 hours.


Figure 5. Toxicological test of chemicals using CCK-8.

Cell line : HeLa
Culture medium : MEM, 10% FBS
Chemicals : Mitomycin-C (MMC),
Dodecylsulfate, sodium salt (SDS)
Incubation : 37°C, 5% CO₂, 2 hours
Detection : 450 nm, Reference : 650 nm

Precautions

1. Since the CCK-8 assay is based on the dehydrogenase activity detection in viable cells, conditions or chemicals that affect dehydrogenase activity in viable cells may cause discrepancy between the actual viable cell number and the cell number determined using the CCK-8 assay.
2. WST-8 may react with reducing agents to generate WST-8 formazan. Please check the background O.D. if reducing agents are used in cytotoxicity assays or cell proliferation assays.
3. Be careful not to introduce bubbles to the wells, since they interfere with the O.D. reading.
4. If the sterilization of the CCK-8 solution is necessary, please filter the solution with a 0.2 μ m membrane.
5. The incubation time varies by the type and number of cells in a well. Generally, leukocytes give weak coloration, thus a long incubation time (up to 4 hours) or a large number of cells ($\sim 10^5$ cells/well) may be necessary.
6. Measure and subtract the O.D. at 600 nm or higher from that of sample if there is a high turbidity in the cell suspension.

Frequently Asked Questions

1. **How many cells should there be in a well?**
For adhesive cells, at least 1000 cells are necessary per well (100 μ l medium). For leukocytes, at least 2500 cells are necessary per well (100 μ l medium) because of low sensitivity. The recommended maximum number of cells per well for the 96-well plate is 25000. If a 24-well or 6-well plate is used for this assay, please calculate the number of cells per well accordingly, and adjust the volume of the CCK-8 solution in a well to 10% of the total volume.
2. **Does CCK-8 stain viable cells?**
No. Since WST-8 and its formazan dye are highly water-soluble, CCK-8 cannot be utilized for cell staining purpose.
3. **Does phenol red affect the assay?**
No. The absorption value of phenol red in a culture medium can be removed by subtracting the absorption value of a blank solution from the absorption value of each well. Therefore, a medium containing phenol red is usable for the CCK-8 assay.
4. **Is CCK-8 toxic to cells?**
Since the toxicity of CCK-8 is so low, the same cells can be used for other cell proliferation assays such as the crystal violet assay, neutral red assay or DNA fluorometric assay after the CCK-8 assay is completed.
5. **I do not have a 450 nm filter. What other filters can I use?**
You can use filters with the absorbance between 430 and 490 nm, even though 450 nm filter gives the best sensitivity.

References

1. M. Ishiyama, Y. Miyazono, K. Sasamoto, Y. Ohkura and K. Ueno, *Talanta*, **1997**, *44*, 1299.
2. H. Tominaga, M. Ishiyama, F. Ohseto, K. Sasamoto, T. Hamamoto, K. Suzuki and M. Watanabe, *Anal. Commun.*, **1999**, *36*, 47.

If you need more information, please contact Dojindo technical service.

Dojindo Laboratories
2025-5 Tabaru, Mashiki-machi, Kamimashiki-gun, Kumamoto
861-2202, Japan Phone: +81-96-286-1515 Fax: +81-96-286-1525
E-mail: info@dojindo.co.jp Web: www.dojindo.co.jp

Dojindo Molecular Technologies, Inc.
Tel: +1-301-987-2667 Web: <http://www.dojindo.com/>
Dojindo EU GmbH
Tel: +49-89-3540-4805 Web: <http://www.dojindo.eu.com/>
Dojindo China Co., Ltd
Tel: +86-21-6427-2302 Web: <http://www.dojindo.cn/>

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CCK-8之Science篇:

1. X-ray screening identifies active site and allosteric inhibitors of SARS-CoV-2 main protease,
Science, **2021**, 372, 642-646
2. Structure-based design of antiviral drug candidates targeting the SARS-CoV-2 main protease,
Science, **2020**, 368, 1331-1335
3. De novo design of potent and resilient hACE2 decoys to neutralize SARS-CoV-2,
Science, **2020**, 370, 1208-1214

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1. Chemical modification of PS-ASO therapeutics reduces cellular protein-binding and improves the therapeutic index, *Nature Biotechnology*, **2019**, 37640-650
2. Multi-omic measurements of heterogeneity in HeLa cells across laboratories, *Nature Biotechnology*, **2019**, 37, 314-322
3. A perfusable, multifunctional epicardial device improves cardiac function and tissue repair, *Nature Medicine*, **2021**, 27, 480-490
4. Dual effects of carbon monoxide on pericytes and neurogenesis in traumatic brain injury, *Nature Medicine*, **2016**, 22, 1335-1341
5. Human induced pluripotent stem cell-derived cardiomyocytes recapitulate the predilection of breast cancer patients to doxorubicin-induced cardiotoxicity, *Nature Medicine*, **2016**, 22(5), 547-556
6. Orphan nuclear receptor NR4A1 regulates transforming growth factor- β signaling and fibrosis, *Nature Medicine*, **2015**, 21(2), 150-158
7. Randomized dose-finding clinical trial of oncolytic immunotherapeutic vaccinia JX-594 in liver cancer, *Nature Medicine*, **2013**, 19, 329-336
8. Compromised CDK1 activity sensitizes BRCA-proficient cancers to PARP inhibition, *Nature Medicine*, **2011**, 17(7), 875-882
9. Mesenchymal stem cell-based tissue regeneration is governed by recipient T lymphocytes via IFN- γ and TNF- α , *Nature Medicine*, **2011**, 17(12), 1594-1601
10. Monoclonal antibody targeting of N-cadherin inhibits prostate cancer growth, metastasis and castration resistance, *Nature Medicine*, **2010**, 16(12), 1414-1420
11. Identification of tendon stem/progenitor cells and the role of the extracellular matrix in their niche, *Nature Medicine*, **2007**, 13(10), 1219-1227
12. Altered TMPRSS2 usage by SARS-CoV-2 Omicron impacts infectivity and fusogenicity, *Nature*, **2022**, doi.org/10.1038/s41586-022-04474-x
13. DAXX represents a new type of protein-folding enabler, *Nature*, **2021**, 597, 132-137
14. DHODH-mediated ferroptosis defence is a targetable vulnerability in cancer, *Nature*, **2021**, 593, 586-590
15. Dietary fructose improves intestinal cell survival and nutrient absorption, *Nature*, **2021**, 597, 263-267
16. Epigenetic therapy inhibits metastases by disrupting premetastatic niches, *Nature*, **2020**, doi.org/10.1038/s41586-020-2054-x
17. ILC2s amplify PD-1 blockade by activating tissue-specific cancer immunity, *Nature*, **2020**, 579, 130-135
18. Potential circadian effects on translational failure for neuroprotection, *Nature*, **2020**, 582, 395-398
19. Chimeric peptidomimetic antibiotics against Gram-negative bacteria, *Nature*, **2019**, doi:10.1038/s41586-019-1665-6
20. PU.1 controls fibroblast polarization and tissue fibrosis, *Nature*, **2019**, 566, 344-349
21. Mitochondrial unfolded protein response controls matrix pre-RNA processing and translation, *Nature*, **2016**, 534, 710-713
22. Transfer of mitochondria from astrocytes to neurons after stroke, *Nature*, **2016**, 535, 551-555
23. Signalling thresholds and negative B-cell selection in acute lymphoblastic leukaemia, *Nature*, **2015**, 521, 357-361
24. The sonic hedgehog factor GLI1 imparts drug resistance through inducible glucuronidation, *Nature*, **2014**, 511(7507), 90-93
25. Sema3A regulates bone-mass accrual through sensory innervations, *Nature*, **2013**, 497(7450), 490-493
26. Acetylation-dependent regulation of endothelial Notch signalling by the SIRT1 deacetylase, *Nature*, **2011**, 473(7346), 234-238
27. Embryonic lethal phenotype reveals a function of TDG in maintaining epigenetic stability, *Nature*, **2011**, 470(7334), 419-423
28. Intravenous delivery of a multi-mechanistic cancer-targeted oncolytic poxvirus in humans, *Nature*, **2011**, 477(7362), 99-102
29. Frequent inactivation of A20 in B-cell lymphomas, *Nature*, **2009**, 459(7247), 712-716
30. Modulation of microRNA processing by p53, *Nature*, **2009**, 460(7254), 529-533
31. Oncogenic mutations of ALK kinase in neuroblastoma, *Nature*, **2008**, 455(7215), 971-974
32. Hyaluronic acid-bilirubin nanomedicine for targeted modulation of dysregulated intestinal barrier, microbiome and immune responses in colitis, *Nature Materials*, **2019**, doi: 10.1038/s41563-019-0462-9
33. Parenchymal and stromal tissue regeneration of tooth organ by pivotal signals reinstated in decellularized matrix, *Nature Materials*, **2019**, 18(6), 627-637
34. Guiding intracortical brain tumour cells to an extracortical cytotoxic hydrogel using aligned polymeric nanofibres, *Nature Materials*, **2014**, 13(3), 308-316
35. Precise targeting of POLR2A as a therapeutic strategy for human triple negative breast cancer, *Nature Nanotechnology*, **2019**, 14, 388-397
36. Metal nanoparticles in the presence of lipopolysaccharides trigger the onset of metal allergy in mice, *Nature Nanotechnology*, **2016**, 11, 808-816

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37. Cross-species chromatin interactions drive transcriptional rewiring in Epstein–Barr virus–positive gastric adenocarcinoma, *Nature Genetics*, **2020**, doi: 10.1038/s41588-020-0665-7
38. Transposable elements drive widespread expression of oncogenes in human cancers, *Nature Genetics*, **2019**, *51*, 611-617
39. Exome-wide analyses identify low-frequency variant in CYP26B1 and additional coding variants associated with esophageal squamous cell carcinoma, *Nature Genetics*, **2018**, *50*, 338-343
40. Pancreatic cancer risk variant in LINC00673 creates a miR-1231 binding site and interferes with PTPN11 degradation, *Nature Genetics*, **2016**, *48*, 747-757
41. Whole-genome mutational landscape and characterization of noncoding and structural mutations in liver cancer, *Nature Genetics*, **2016**, *48*, 500-509
42. A recurrent inactivating mutation in RHOA GTPase in angioimmunoblastic T cell lymphoma, *Nature Genetics*, **2014**, *46*, 371-375
43. Exome sequencing identifies somatic gain-of-function PPM1D mutations in brainstem gliomas, *Nature Genetics*, **2014**, *46*, 726-730
44. Somatic RHOA mutation in angioimmunoblastic T cell lymphoma, *Nature Genetics*, **2014**, *46*, 171-175
45. Genome-wide association study identifies common variants in SLC39A6 associated with length of survival in esophageal squamous-cell carcinoma, *Nature Genetics*, **2013**, *45*(6), 632-638
46. Recurrent mutations in multiple components of the cohesin complex in myeloid neoplasms, *Nature Genetics*, **2014**, *45*, 1232-1237
47. Common variants at 11q12, 10q26 and 3p11.2 are associated with prostate cancer susceptibility in Japanese, *Nature Genetics*, **2012**, *44*(4), 426-429
48. Exome sequencing of hepatitis B virus-associated hepatocellular carcinoma, *Nature Genetics*, **2012**, *44*(10), 1117-1121
49. LIN28B induces neuroblastoma and enhances MYCN levels via let-7 suppression, *Nature Genetics*, **2012**, *44*(11), 1199-1206
50. Heterozygosity with respect to Zfp148 causes complete loss of fetal germ cells during mouse embryogenesis, *Nature Genetics*, **2003**, *33*, 172-176
51. The tumor suppressor PTEN has a critical role in antiviral innate immunity, *Nature Immunology*, **2016**, *17*, 241-249
52. Neddylation of PTEN regulates its nuclear import and promotes tumor development, *Cell Research*, **2021**, *31*, 291-311
53. Pericytes augment glioblastoma cell resistance to temozolomide through CCL5-CCR5 paracrine signaling, *Cell Research*, **2021**, <https://doi.org/10.1038/s41422-021-00528-3>
54. SIDT1-dependent absorption in the stomach mediates host uptake of dietary and orally administered microRNAs, *Cell Research*, **2021**, *31*, 247-258
55. Inhibition of SARS-CoV-2 (previously 2019-nCoV) infection by a highly potent pan-coronavirus fusion inhibitor targeting its spike protein that harbors a high capacity to mediate membrane fusion, *Cell Research*, **2020**, *30*, 343-355
56. The role of ferroptosis in ionizing radiation-induced cell death and tumor suppression, *Cell Research*, **2020**, *30*, 146-162
57. A two-step lineage reprogramming strategy to generate functionally competent human hepatocytes from fibroblasts, *Cell Research*, **2019**, *29*, 696-710
58. Expansion and differentiation of human hepatocyte-derived liver progenitor-like cells and their use for the study of hepatotropic pathogens, *Cell Research*, **2019**, *29*, 8-22
59. PARylation regulates stress granule dynamics, phase separation, and neurotoxicity of disease-related RNA-binding proteins, *Cell Research*, **2019**, *29*, 233-247
60. A dynamic N6-methyladenosine methylome regulates intrinsic and acquired resistance to tyrosine kinase inhibitors, *Cell Research*, **2018**, *28*, 1062-1076
61. An NF90/NF110-mediated feedback amplification loop regulates dicer expression and controls ovarian carcinoma progression, *Cell Research*, **2018**, *28*, 556-571
62. Integrated genomic analysis identifies deregulated JAK/STAT-MYC-biosynthesis axis in aggressive NK-cell leukemia, *Cell Research*, **2018**, *28*, 172-186
63. Existing drugs as broad-spectrum and potent inhibitors for Zika virus by targeting NS2B-NS3 interaction, *Cell Research*, **2017**, *27*, 1046-1064
64. Epigenetic silencing of microRNA-149 in cancer-associated fibroblasts mediates prostaglandin E2/interleukin-6 signaling in the tumor microenvironment, *Cell Research*, **2015**, *25*, 588-603

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1. Virological characteristics of the SARS-CoV-2 Omicron BA.2 spike, *Cell*, **2022**, doi.org/10.1016/j.cell.2022.04.035
2. In Vivo CRISPR Screens Identify E3 Ligase Cop1 as a Modulator of Macrophage Infiltration and Cancer Immunotherapy Target, *Cell*, **2021**, DOI:https://doi.org/10.1016/j.cell.2021.09.006
3. A Bacterial Effector Reveals the V-ATPase-ATG16L1 Axis that Initiates Xenophagy, *Cell*, **2019**, 178(3), 552-566
4. Higher-Order Clustering of the Transmembrane Anchor of DR5 Drives Signaling, *Cell*, **2019**, 176(6), 1477-1489
5. Insulin Receptor Associates with Promoters Genome-wide and Regulates Gene Expression, *Cell*, **2019**, 177(3), 722-736
6. B-Cell-Specific Diversion of Glucose Carbon Utilization Reveals a Unique Vulnerability in B Cell Malignancies, *Cell*, **2018**, 173(2), 470-484
7. *Fusobacterium nucleatum* Promotes Chemoresistance to Colorectal Cancer by Modulating Autophagy, *Cell*, **2017**, 170(3), 548-563
8. Methyltransferase SETD2-Mediated Methylation of STAT1 Is Critical for Interferon Antiviral Activity, *Cell*, **2017**, 170(3), 492-506
9. Loss of 5-Hydroxymethylcytosine Is an Epigenetic Hallmark of Melanoma, *Cell*, **2012**, 150(6), 1135-1146
10. Mapping the Hallmarks of Lung Adenocarcinoma with Massively Parallel Sequencing, *Cell*, **2012**, 150(6), 1107-1120
11. Iron-Export Ferroxidase Activity of β -Amyloid Precursor Protein Is Inhibited by Zinc in Alzheimer's Disease, *Cell*, **2010**, 142(6), 857-867
12. Autocrine VEGF Signaling Is Required for Vascular Homeostasis, *Cell*, **2007**, 130(4), 691-703
13. Genomic and Transcriptomic Characterization of Natural Killer T Cell Lymphoma, *Cancer Cell*, **2020**, 37, 403-419
14. CDK7 Inhibition Potentiates Genome Instability Triggering Anti-tumor Immunity in Small Cell Lung Cancer, *Cancer Cell*, **2020**, 37, 37-54
15. Wild-Type p53 Promotes Cancer Metabolic Switch by Inducing PUMA-Dependent Suppression of Oxidative Phosphorylation, *Cancer Cell*, **2019**, 35(2), 191-203
16. Genomic and Epigenomic Profiling of High-Risk Intestinal Metaplasia Reveals Molecular Determinants of Progression to Gastric Cancer, *Cancer Cell*, **2018**, 33(1), 137-150
17. Transcriptional Regulation of the Warburg Effect in Cancer by SIX1, *Cancer Cell*, **2018**, 33(3), 368-385
18. Enhanced Fructose Utilization Mediated by SLC2A5 Is a Unique Metabolic Feature of Acute Myeloid Leukemia with Therapeutic Potential, *Cancer Cell*, **2016**, 30(5), 779-791
19. Targeting Transcriptional Addictions In Small Cell Lung Cancer With a Covalent CDK7 Inhibitor, *Cancer Cell*, **2014**, 26(6), 909-922
20. Nrf2 Redirects Glucose and Glutamine into Anabolic Pathways in Metabolic Reprogramming, *Cancer Cell*, **2012**, 22(1), 66-79
21. Polycomb-Mediated Loss of miR-31 Activates NIK-Dependent NF- κ B Pathway in Adult T Cell Leukemia and Other Cancers, *Cancer Cell*, **2012**, 21(1), 121-135
22. E2F1-Regulated MicroRNAs Impair TGF β -Dependent Cell-Cycle Arrest and Apoptosis in Gastric Cancer, *Cancer Cell*, **2008**, 13, 272-286
23. An Upstream Open Reading Frame in Phosphatase and Tensin Homolog Encodes a Circuit Breaker of Lactate Metabolism, *Cell Metabolism*, **2021**, 33, 28-144
24. Apaf-1 Pyroptosome Senses Mitochondrial Permeability Transition, *Cell Metabolism*, **2021**, 33, 1-13
25. Metabolic-Pathway-Based Subtyping of Triple[1]Negative Breast Cancer Reveals Potential Therapeutic Targets, *Cell Metabolism*, **2021**, 33, 1-14
26. A Novel Allosteric Inhibitor of Phosphoglycerate Mutase 1 Suppresses Growth and Metastasis of Non-Small-Cell Lung Cancer, *Cell Metabolism*, **2019**, 30(6), 1107-1119
27. IGF-2 Preprograms Maturing Macrophages to Acquire Oxidative Phosphorylation-Dependent Anti-inflammatory Properties, *Cell Metabolism*, **2019**, 29(6), 1363-1375
28. Liver damage, inflammation, and enhanced tumorigenesis after persistent mTORC1 inhibition, *Cell Metabolism*, **2014**, 20(1), 133-44
29. Loss of SATB1 Induces p21-Dependent Cellular Senescence in Post-mitotic Dopaminergic Neurons, *Cell Stem Cell*, **2019**, 25(4), 514-530
30. PPM1D Mutations Drive Clonal Hematopoiesis in Response to Cytotoxic Chemotherapy, *Cell Stem Cell*, **2018**, 23(5), 700-713
31. Lipid Desaturation Is a Metabolic Marker and Therapeutic Target of Ovarian Cancer Stem Cells, *Cell Stem Cell*, **2017**, 20(3), 303-314

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1. The effect of novel acridine-based agents with topoisomerase II inhibitor on mesothelioma cell proliferation and apoptosis, *Journal of Clinical Oncology*, **2011**, 15, DOI:10.1200/jco.2011.29.15_suppl.e13507
2. Genetic risk of extranodal natural killer T-cell lymphoma: a genome-wide association study in multiple populations, *Lancet Oncology*, **2019**, pii: S1470-2045(19)30799-5
3. A Secondary Mutation in BRAF Confers Resistance to RAF Inhibition in a BRAF V600E-Mutant Brain Tumor, *Cancer Discovery*, **2018**, 8(9), 1130-1141
4. Epigenomic Promoter Alterations Amplify Gene Isoform and Immunogenic Diversity in Gastric Adenocarcinoma, *Cancer Discovery*, **2017**, 7(6), 630-651
5. A novel photothermally controlled multifunctional scaffold for clinical treatment of osteosarcoma and tissue regeneration, *Materials Today*, **2020**, doi.org/10.1016/j.mattod.2019.12.005
6. 3D Heterogeneous Device Arrays for Multiplexed Sensing Platforms Using Transfer of Perovskites, *Advanced Materials*, **2021**, 33, 2101093
7. Hydroxyapatite Nanorods Function as Safe and Effective Growth Factors Regulating Neural Differentiation and Neuron Development, *Advanced Materials*, **2021**, 2100895
8. Self-Activatable Photo-Extracellular Vesicle for Synergistic Trimodal Anticancer Therapy, *Advanced Materials*, **2021**, 33, 2005562
9. Stretchable and Highly Permeable Nanofibrous Sensors for Detecting Complex Human Body Motion, *Advanced Materials*, **2021**, 2102488
10. A Graphdiyne Oxide-Based Iron Sponge with Photothermally Enhanced Tumor-Specific Fenton Chemistry, *Advanced Materials*, **2020**, 2000038
11. A Highly Efficient Tumor-Targeting Nanoprobe with a Novel Cell Membrane Permeability Mechanism, *Advanced Materials*, **2019**, 31(12), e1807456
12. A Self-Pumping Dressing for Draining Excessive Biofluid around Wounds, *Advanced Materials*, **2019**, 31(5), e1804187
13. Biomimetic Metal-Organic Framework Nanoparticles for Cooperative Combination of Antiangiogenesis and Photodynamic Therapy for Enhanced Efficacy, *Advanced Materials*, **2019**, 31(15), 1808200
14. CD44-Specific A6 Short Peptide Boosts Targetability and Anticancer Efficacy of Polymersomal Epirubicin to Orthotopic Human Multiple Myeloma, *Advanced Materials*, **2019**, doi.org/10.1002/adma.201904742
15. Chirality Controls Mesenchymal Stem Cell Lineage Diversification through Mechanoresponses, *Advanced Materials*, **2019**, doi.org/10.1002/adma.201900582
16. Engineering Biomimetic Plasmosomes for pH-Responsive Drug Delivery and Enhanced Antitumor Activity, *Advanced Materials*, **2019**, 31(32), e1900795
17. Highly Efficient and Environmentally Friendly Fabrication of Robust, Programmable, and Biocompatible Anisotropic, All-Cellulose, Wrinkle-Patterned Hydrogels for Cell Alignment, *Advanced Materials*, **2019**, 31(46), e1904762
18. Securing the Payload, Finding the Cell, and Avoiding the Endosome: Peptide-Targeted, Fusogenic Porous Silicon Nanoparticles for Delivery of siRNA, *Advanced Materials*, **2019**, 31(35), e1902952
19. Semiconducting Polymer Nanoparticles for Centimeters-Deep Photoacoustic Imaging in the Second Near-Infrared Window, *Advanced Materials*, **2017**, 29(41), 1703403
20. Structurally and Functionally Optimized Silk-Fibroin-Gelatin Scaffold Using 3D Printing to Repair Cartilage Injury In Vitro and In Vivo, *Advanced Materials*, **2017**, 29(29), doi: 10.1002/adma.201701089
21. DNA-Nanostructure-Gold-Nanorod Hybrids for Enhanced In Vivo Photoacoustic Imaging and Photothermal Therapy, *Advanced Materials*, **2016**, 28(45):10000-10007

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